


## ASSEMBLY INSTRUCTIONS: WOOD SCIENCE TABLE (WST)



HARDWARE		
QUANTITY	PART NUMBER	DESCRIPTION
4	1010528	PERIPHERAL GLIDE
8	10105661	3/8 X 2" HEX HEAD LAG BOLT

### INSTRUCTIONS

1. Lay the top upside down on a clean, scratch free surface.
2. Place leg in corner.
3. Attach leg to top with two 3/8 X 2" hex head lag bolts per leg.
4. Turn table upright.